



Product Datasheet

MET antibody (orb127420)

Description	Rabbit polyclonal antibody to MET
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P08581
Dilution Range	WB: 1:500-1:3000, ELISA: 1:10000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from K...



Western blot analysis of extracts from J...