

## Product Datasheet

# Recombinant-SARS-CoV-2 (COVID-19) RBD Neutralizing Antibody [S2H97] (orb1274209)

**Description** Recombinant-SARS-CoV-2 (COVID-19)  
RBD Neutralizing Antibody [S2H97]

**Species/Host** Human

**Reactivity** Virus

**Conjugation** Unconjugated

**Tested Applications** ELISA, WB

**Preservatives** PBS,pH7.4

**Form/Appearance** Liquid

**Concentration** batch dependent

**Storage** Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for 1-2 weeks. Store at -20 °C for 12 months. Store at -80 °C for long term storage.

**Note** For research use only


**Application notes** WB:1:1000-2000; ELISA:1:5000-1:10000

**Isotype** IgG1-kappa

**Clonality** Recombinant

**Dilution Range** WB:1:1000-2000; ELISA:1:5000-1:10000

**Expiration Date** 12 months from date of receipt.



SDS-  
PAGE  
image