

Product Datasheet

Recombinant-SARS-CoV-2 (COVID-19) RBD Neutralizing Antibody [S2E12] (orb1274210)

Description Recombinant-SARS-CoV-2 (COVID-19)
RBD Neutralizing Antibody [S2E12]

Species/Host Human

Reactivity Virus

Conjugation Unconjugated

Tested Applications ELISA, WB

Preservatives PBS,pH7.4

Form/Appearance Liquid

Concentration batch dependent

Storage Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for 1-2 weeks. Store at -20 °C for 12 months. Store at -80 °C for long term storage.

Note For research use only


Application notes WB:1:1000-2000; ELISA:1:5000-1:10000

Isotype IgG1-kappa

Clonality Recombinant

Dilution Range WB:1:1000-2000; ELISA:1:5000-1:10000

Expiration Date 12 months from date of receipt.



SDS-
PAGE
image