

Product Datasheet

SARS-CoV-2 Nucleocapsid Antibody, monoclonal (orb1274212)

Description

SARS-CoV-2 Nucleocapsid Antibody, monoclonal

Reactivity

Human

Conjugation

Unconjugated

Tested

Applications

ELISA, WB

Immunogen

Anti-SARS-CoV-2 Nucleocapsid Antibody, Human IgG1 is isolated from a SARS-CoV-2 infected patient and is recombinantly produced from human 293 cells (HEK293).

Target

N

Preservatives

Lyophilized from 0.22 μ m filtered solution in PBS, pH 7.4 . Normally trehalose is added as protectant before lyophilization.

Form/Appearance

Lyophilized

Concentration

batch dependent

Storage

For long term storage, the product should be stored at lyophilized state at -20°C or lower..Please avoid repeated freeze-thaw cycles.This product is stable after storage at:-20 to -70°C for 12 months in lyophilized state from date of receipt;.-70°C for 3 months under sterile conditions after reconstitution.

Note

For research use only

Application notes

This antibody can be paired with other Anti-SARS-CoV-2 Nucleocapsid antibodies to detect SARS-CoV-2 Nucleocapsid protein in sandwich ELISA or LFA assay.

Isotype

Human IgG1/kappa

Clonality

Recombinant

Uniprot ID

P0DTC9


Dilution Range

This antibody can be paired with other Anti-SARS-CoV-2 Nucleocapsid antibodies to detect SARS-CoV-2 Nucleocapsid protein in sandwich ELISA or LFA assay.


Expiration Date

12 months from date of receipt.

kDa M R NR



Anti-SARS-CoV-2 Nucleocapsid Antibody, H...



Immobilized SARS-CoV-2 Nucleocapsid prot...