

Product Datasheet

**SARS-CoV-2 Spike NTD Neutralizing Antibody,
monoclonal, Chimeric (orb1274218)**

Description

SARS-CoV-2 Spike NTD Neutralizing Antibody,



Reactivity

Mouse

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Anti-SARS-CoV-2 Spike NTD Antibody, Chimeric mAb is a chimeric monoclonal antibody combining the constant domains of the human IgG1 molecule with mouse variable regions.


Anti-SARS-CoV-2 Spike NTD Antibody, Chim...

Target

S

Preservatives

Lyophilized from 0.22 μ m filtered solution in PBS, pH 7.4 . Normally trehalose is added as protectant before lyophilization.



Form/Appearance

Lyophilized


Concentration

batch dependent

Storage

For long term storage, the product should be stored at lyophilized state at -20°C or lower..Please avoid repeated freeze-thaw cycles.This product is stable after storage at:-20 to -70°C for 12 months in lyophilized state from date of receipt;-70°C for 3 months under sterile conditions after reconstitution.

Immobilized SARS-CoV-2 S1 protein NTD, H...



Note

For research use only

Isotype

Human IgG1/kappa

Clonality

Recombinant

Uniprot ID

P0DTC2

Anti-SARS-CoV-2 Spike NTD Antibody, Chim...

Expiration Date

12 months from date of receipt.