


Product Datasheet

PTH Antibody (orb1274225)

Description	PTH Antibody	
Species/Host	Rabbit	
Reactivity	Human	IHC testing of FFPE human parathyroid gl...
Conjugation	Unconjugated	
Tested Applications	FC, IF, IHC-P	
Immunogen	A synthetic peptide from the N-terminal region of human Parathyroid hormone was used as the immunogen for the recombinant PTH antibody.	
Target	PTH	
Preservatives	PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide	
Form/Appearance	Liquid	
Concentration	0.2 mg/mL	
Storage	Aliquot and Store at 2-8°C. Avoid freez-thaw cycles.	
Note	For research use only	
Application notes	Flow Cytometry: 0.5-1 ug/million cells in 0.1ml Immunofluorescence: 0.5-1 ug/ml Immunohistochemistry (FFPE): 0.5-1 ug/ml for 30 min at RT Prediluted IHC only format: incubate for 30 min at RT (1)Optimal dilution of the recombinant PTH antibody should be determined by the researcher.1. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and incubate at RT for 30 min.	
Isotype	IgG, kappa	
Clonality	Recombinant	
Uniprot ID	P01270	
Dilution Range	Flow Cytometry: 0.5-1 ug/million cells in 0.1ml Immunofluorescence: 0.5-1 ug/ml Immunohistochemistry (FFPE): 0.5-1 ug/ml for 30 min at RT Prediluted IHC only format: incubate for 30 min at RT (1)Optimal dilution of the recombinant PTH antibody should be determined by the researcher.1. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and incubate at RT for 30 min.	
Expiration Date	12 months from date of receipt.	