

Product Datasheet

MLANA Antibody (orb1274227)

Description

MLANA Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

FC, IF, IHC-P, WB

Applications

Immunogen

Human protein was used as the immunogen for this recombinant Melan-A antibody.

Target

MLANA

Preservatives

PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide

Form/Appearance

Liquid

Concentration

0.2 mg/mL

Storage

Aliquot and Store at 2-8°C. Avoid freeze-thaw cycles.

Note

For research use only

Application notes

Western blot: 0.5-1 ug/ml
Flow Cytometry: 0.5-1ug/10⁶ cells
Immunofluorescence: 0.5-1 ug/ml
Immunohistochemistry (FFPE): 0.5-1 ug/ml for 30 minutes at RT
Prediluted IHC only format: incubate for 30 min at RT (1)
The optimal dilution of the recombinant Melan-A antibody for each application should be determined by the researcher.
1. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and incubate at RT for 30 min.

Isotype

IgG, k

Clonality

Recombinant

Uniprot ID


[Q16655](#)

Dilution Range


Western blot: 0.5-1 ug/ml
Flow Cytometry: 0.5-1ug/10⁶ cells
Immunofluorescence: 0.5-1 ug/ml
Immunohistochemistry (FFPE): 0.5-1 ug/ml for 30 minutes at RT
Prediluted IHC only format: incubate for 30 min at RT (1)
The optimal dilution of the recombinant Melan-A antibody for each application should be determined by the researcher.
1. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and incubate at RT for 30 min.

Expiration Date

12 months from date of receipt.



IHC testing
FFPE human
skin with
recombi...



IHC testing of
FFPE human
melanoma
with ...