

Product Datasheet

TNFRSF8 Antibody (orb1274238)

Description

TNFRSF8 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IF, IHC-P

Applications

Immunogen

A partial human protein was used as the immunogen for the recombinant CD30 antibody.

Target

TNFRSF8

Preservatives

PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide

Form/Appearance

Liquid

Concentration

0.2 mg/mL

Storage

Aliquot and Store at 2-8°C. Avoid freez-thaw cycles.

Note

For research use only

Application notes

Flow Cytometry: 0.5-1 ug/million cells in 0.1ml
Immunofluorescence: 1-2 ug/ml
Immunohistochemistry (FFPE): 0.5-1 ug/ml for 30 min at RT
Prediluted IHC only format: incubate for 30 min at RT
(1)Optimal dilution of the recombinant CD30 antibody should be determined by the researcher.1. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and incubate at RT for 30 min.

Isotype

IgG, kappa

Clonality

Recombinant

Uniprot ID


P28908

Dilution Range

Flow Cytometry: 0.5-1 ug/million cells in 0.1ml
Immunofluorescence: 1-2 ug/ml
Immunohistochemistry (FFPE): 0.5-1 ug/ml for 30 min at RT
Prediluted IHC only format: incubate for 30 min at RT
(1)Optimal dilution of the recombinant CD30 antibody should be determined by the researcher.1. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and incubate at RT for 30 min.

Expiration Date

12 months from date of receipt.



IHC testing of FFPE human Hodgkin's lymphoma tissue.