

Product Datasheet

MS4A1 Antibody (orb1274246)

Description

MS4A1 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC

Applications

Immunogen

Full length human CD20 protein was used as the immunogen for this recombinant CD20 antibody.

Target

MS4A1

Preservatives

PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide

Form/Appearance

Liquid

Concentration

0.2 mg/mL

Storage

Aliquot and Store at 2-8°C. Avoid freeze-thaw cycles.

Note

For research use only

Application notes

Immunohistochemistry (FFPE): 0.5-1 ug/ml for 30 minutes at RT
 Prediluted IHC only format: incubate for 30 min at RT
 (1) Titration of the recombinant CD20 antibody may be required for optimal performance.1. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and incubate at RT for 30 min.

Isotype

IgG, kappa

Clonality

Recombinant

Uniprot ID


P11836

Dilution Range

Immunohistochemistry (FFPE): 0.5-1 ug/ml for 30 minutes at RT
 Prediluted IHC only format: incubate for 30 min at RT
 (1) Titration of the recombinant CD20 antibody may be required for optimal performance.1. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and incubate at RT for 30 min.

Expiration Date

12 months from date of receipt.



IHC testing of recombinant CD20 antibody...