

Product Datasheet

PECAM1 Antibody (orb1274252)

Description

PECAM1 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IF, IHC

Applications

Immunogen

A full length recombinant human protein was used as the immunogen for the rabbit monoclonal CD31 antibody.

Target

PECAM1

Preservatives

PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide

Form/Appearance

Liquid

Concentration

0.2 mg/mL

Storage

Aliquot and Store at 2-8°C. Avoid freez-thaw cycles.

Note

For research use only

Application notes

Immunofluorescence: 0.5-1 ug/ml Immunohistochemistry (FFPE): 0.5-1 ug/ml for 30 minutes at RT (1)Titering of the CD31 antibody may be required for optimal performance.1. Staining of formalin-fixed tissues is enhance by boiling sections in 1mM EDTA, pH 7.5-8.5, for 10-20 min followed by cooling at RT for 20 minutes.

Isotype

IgG, k

Clonality

Recombinant

Uniprot ID


P16284

Dilution Range

Immunofluorescence: 0.5-1 ug/ml Immunohistochemistry (FFPE): 0.5-1 ug/ml for 30 minutes at RT (1)Titering of the CD31 antibody may be required for optimal performance.1. Staining of formalin-fixed tissues is enhance by boiling sections in 1mM EDTA, pH 7.5-8.5, for 10-20 min followed by cooling at RT for 20 minutes.

Expiration Date

12 months from date of receipt.



IHC staining of
FFPE human
angiocarcinom...