

Product Datasheet

WT1 Antibody (orb1274254)

Description

WT1 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IF, IHC-P

Applications

Immunogen

A recombinant protein was used as the immunogen for the WT1 antibody.

Target

WT1

Preservatives

PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide

Form/Appearance

Liquid

Concentration

0.2 mg/mL

Storage

Aliquot and Store at 2-8°C. Avoid freeze-thaw cycles.

Note

For research use only

Application notes

Immunofluorescence: 0.5-1 ug/ml
Immunohistochemistry (FFPE): 0.5-1 ug/ml for 30 min at RT (1)
Optimal dilution of the WT1 antibody should be determined by the researcher.
1. Staining of formalin-fixed tissues requires boiling tissue sections in 10mM Citrate buffer, pH 6.0, for 10-20 min followed by cooling at RT for 20 min.

Isotype

IgG, k

Clonality

Recombinant

Uniprot ID


[P19544](#)

Dilution Range

Immunofluorescence: 0.5-1 ug/ml
Immunohistochemistry (FFPE): 0.5-1 ug/ml for 30 min at RT (1)
Optimal dilution of the WT1 antibody should be determined by the researcher.
1. Staining of formalin-fixed tissues requires boiling tissue sections in 10mM Citrate buffer, pH 6.0, for 10-20 min followed by cooling at RT for 20 min.

Expiration Date

12 months from date of receipt.



Formalin-fixed,
paraffin-embedded
human ...