

Product Datasheet

CFTR Antibody (orb1274272)

Description

CFTR Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IF, IHC-P

Applications

Immunogen

A recombinant partial protein corresponding to amino acids 258-385 from the human protein was used as the immunogen for this recombinant CFTR antibody.

Target

CFTR

Preservatives

PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide

Form/Appearance

Liquid

Concentration

0.2 mg/mL

Storage

Aliquot and Store at 2-8°C. Avoid freeze-thaw cycles.

Note

For research use only

Application notes

IF: 0.5-1 ug/ml IHC (FFPE): 1-2 ug/ml for 30 min at RT
Prediluted IHC only format: incubate for 30 min at RT
(1) The concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the recombinant CFTR antibody to be titrated up or down for optimal performance. 1. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and incubate at RT for 30 min.

Isotype

IgG, k

Clonality

Recombinant

Uniprot ID


[P13569](#)

Dilution Range

IF: 0.5-1 ug/ml IHC (FFPE): 1-2 ug/ml for 30 min at RT
Prediluted IHC only format: incubate for 30 min at RT
(1) The concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the recombinant CFTR antibody to be titrated up or down for optimal performance. 1. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and incubate at RT for 30 min.

Expiration Date

12 months from date of receipt.



IHC testing of
FFPE human
pancreas
with ...