


Product Datasheet

SARS-CoV-2 Spike S1 Control Antibody [4D1J6], Chimeric [Out of stock] (orb1274384)

Description	SARS-CoV-2 Spike S1 Control Antibody [4D1J6], Chimeric [Out of stock]
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	HEK293 derived SARS-COV-2 SPIKE S1, amino acids Val16-Arg685
Target	S
Preservatives	Supplied as a 0.2um filtered solution in PBS, pH 7.4.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	This antibody can be stored at 2°C to -8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Chimeric MAb
Clonality	Recombinant
Uniprot ID	P0DTC2
NCBI	YP_009724390.1
Expiration Date	12 months from date of receipt.



This standard curve is only for demonstr...