


Product Datasheet

**SARS-COV-2 Nucleocapsid Protein Antibody
[9A4B3] [Out of stock] (orb1274386)**


Description

SARS-COV-2 Nucleocapsid Protein Antibody


Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	E.coli derived SARS-COV-2 Nucleoproteins Ser2-Ala419-Gly335Ala Accession #QHD43423.2
Target	N
Preservatives	Lyophilized from a 0.2 um filtered solution in PBS, pH 7.4 with antibody at a final concentration of 1 mg/ml.
Form/Appearance	Lyophilized
Concentration	1 mg/ml
Storage	2°C - 8°C for one month, -20°C to -80°C for one year. Avoid repeated freeze-thaw cycles
Note	For research use only
Application notes	SARS-COV-2 Nucleocapsid Sandwich Elisa: ELISA Capture: Recommended Concentration:0.5-2ug/mL, Sample: SARS-CoV-2 Nucleocapsid Protein Antibody-CAP) ELISA Capture: Recommended Concentration:0.005-0.02ug/mL, Sample: SARS-CoV-2 Nucleocapsid Protein Antibody Standard:Recommended Concentration: 156-10000pg/mL, Sample: Recombinant SARS-COV-2 Nucleocapsid Protein with His tag
Isotype	Rabbit IgG
Clonality	Recombinant
Uniprot ID	P0DTC9
NCBI	QHD43423.2
Dilution Range	SARS-COV-2 Nucleocapsid Sandwich Elisa: ELISA Capture: Recommended Concentration:0.5-2ug/mL, Sample: SARS-CoV-2 Nucleocapsid Protein Antibody-CAP) ELISA Capture: Recommended Concentration:0.005-0.02ug/mL, Sample: SARS-CoV-2 Nucleocapsid Protein Antibody Standard:Recommended Concentration: 156-10000pg/mL, Sample: Recombinant SARS-COV-2 Nucleocapsid Protein



This standard curve is only for demonstr...



No cross-reactivity in ELISA with MERS-C...



Recombinant SARS-CoV-2 Nucleocapsid Prot...