


Product Datasheet

**SARS-COV-2 Nucleocapsid Protein Antibody
(biotin) [1B1A9] [Out of stock] (orb1274394)**


Description

SARS-COV-2 Nucleocapsid Protein Antibody


Species/Host	Rabbit
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Hek293 derived SARS-COV-2 Nucleocapsid Ser2-Ala419 Accession #YP_009724397.2
Target	N
Preservatives	Supplied as a 0.2um filtered solution in PBS, pH 7.4 containing SARS-CoV-2 (2019-nCoV) Nucleocapsid Antibody at a concentration of 1mg/mL.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	This antibody can be stored at 2°C - 8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	SARS-COV-2 Nucleocapsid Sandwich Elisa: ELISA Capture: Recommended Concentration: 1-4ug/mL, Sample: SARS-CoV-2 Nucleocapsid Protein Antibody (CAP) ELISA Detection: Recommended Concentration: 0.05-0.2ug/mL, Sample: Biotinylated SARS-CoV-2 Nucleocapsid Protein Antibody Standard: Recommended Concentration: 1.56-100pg/mL, Sample: Recombinant SARS-COV-2 Nucleocapsid Protein with His tag
Isotype	Rabbit IgG
Clonality	Recombinant
Uniprot ID	P0DTC9
NCBI	YP_009724397.2
Dilution Range	SARS-COV-2 Nucleocapsid Sandwich Elisa: ELISA Capture: Recommended Concentration: 1-4ug/mL, Sample: SARS-CoV-2 Nucleocapsid Protein Antibody (CAP) ELISA Detection:



This standard curve is only for demonstr...



No cross-reactivity in ELISA with MERS-C...



Cat. No. orb1274386, captured on COOH ch...