


Product Datasheet

AKT2 Antibody (orb1274477)

Description	AKT2 Antibody
Species/Host	Mouse
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ChIP, ICC, IF, IHC, IP, WB
Immunogen	Recombinant protein of human AKT2
Target	AKT2
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:500 - 1:2000IHC: 1:50 - 1:100ICC: 1:50 - 1:200IF: 1:50 - 1:100IP: 1:20 - 1:50ChIP: 1:20 - 1:100
Isotype	IgG
Clonality	Monoclonal
MW	Observed: 60kDa
Uniprot ID	P31751
Dilution Range	WB: 1:500 - 1:2000IHC: 1:50 - 1:100ICC: 1:50 - 1:200IF: 1:50 - 1:100IP: 1:20 - 1:50ChIP: 1:20 - 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from n...