



Product Datasheet

PARP1 Antibody (orb1274565)


Description	PARP1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ChIP, IF, IHC, IP, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 81-390 of human PARP1 (NP_001609.2).
Target	PARP1
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Application notes	WB: 1:1000 - 1:3000IHC: 1:50 - 1:200IF: 1:50 - 1:200IP: 1:50 - 1:200CHIP: 1:50 - 1:200
Isotype	IgG
Clonality	Polyclonal
MW	Observed: 100kDa, 110kDa, 125kDa
Uniprot ID	P09874
Dilution Range	WB: 1:1000 - 1:3000IHC: 1:50 - 1:200IF: 1:50 - 1:200IP: 1:50 - 1:200CHIP: 1:50 - 1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Western blot analysis of extracts of Jur...



Western blot analysis of extracts from n...