



Product Datasheet

Recombinant SARS-CoV-2 (COVID-19) RBD Neutralizing antibody [AbA128] (orb1274586)

Description	Recombinant SARS-CoV-2 (COVID-19) RBD Neutralizing antibody [AbA128]
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant SARS-CoV-2 RBD protein
Preservatives	PBS, pH7.5
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for 1-2 weeks. Store at -20 °C for 12 months. Store at -80 °C for long term storage.
Note	For research use only
Application notes	WB:1:1000-2000ELISA:1:5000-1:10000
Isotype	Human IgG
Clonality	Recombinant
Dilution Range	WB:1:1000-2000ELISA:1:5000-1:10000
Expiration Date	12 months from date of receipt.



Anti-RBD
antibody
(orb1274586)
was used ...



Coomassie
blue staining
non-reduced
and ...