



Product Datasheet

Recombinant SARS-CoV-2 (COVID-19) (S1) -1 Antibody [CR3022] (orb1274599)

Description	Recombinant SARS-CoV-2 (COVID-19) (S1) -1 Antibody [CR3022]
Species/Host	Human
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant SARS-CoV-2 S1 protein
Preservatives	PBS,pH7.5
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for 1-2 weeks. Store at -20 °C for 12 months. Store at -80 °C for long term storage.
Note	For research use only
Application notes	ELISA:1:5000-1:10000
Isotype	Human IgG1
Clonality	Recombinant
Dilution Range	ELISA:1:5000-1:10000
Expiration Date	12 months from date of receipt.



Coomassie
blue
staining
non-
reduced
and ...



Sandwich
ELISA for
SARS-CoV-
2 (COVID-
19)...