

Product Datasheet

SARS-CoV-2 (COVID-19) Nucleocapsid Antibody [85C10] (orb1274604)

Description SARS-CoV-2 (COVID-19) Nucleocapsid Antibody [85C10]

Species/Host Rabbit

Reactivity Virus

Conjugation Unconjugated

Tested Applications ELISA, WB

Immunogen SARS-CoV-2 nucleoprotein (NP) specific peptide

Target N

Preservatives 1X PBS, 0.02% NaN₃

Form/Appearance Liquid

Concentration batch dependent

Storage SARS-CoV-2 (COVID-19) Nucleoprotein Antibody [85C10] can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note For research use only

Application notes ELISA


Isotype IgG1k

Clonality Recombinant


NCBI [QHD43423](#)

Dilution Range ELISA

Expiration Date 12 months from date of receipt.



Microtiter wells were coated with SARS-CoV-2...



Western blot analysis of SARS-CoV-2 and S...