


Product Datasheet

SARS-CoV-2 (COVID-19) Nucleocapsid Antibody [84D7] (orb1274606)


Description

SARS-CoV-2 (COVID-19) Nucleocapsid


COVID-19 (+) COVID-19 (-)



Immunohistochemistry
Validation of SARS-...



Microtiter wells were
coated with SARS-C...



Western blot analysis
of SARS-CoV2 and S...

Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	SARS-CoV-2 nucleoprotein (NP) specific peptide
Target	N
Preservatives	1X PBS, 0.02% NaN ₃ , 50% Glycerol, 0.1% rAlbumin
Form/Appearance	Liquid
Concentration	batch dependent
Storage	SARS-CoV-2 (COVID-19) Nucleoprotein Antibody [84D7] can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	Paired with 75G5a, used for specific, high-sensitivity detection of SARS-CoV-2 nucleocapsid protein (NP).
Isotype	IgG1k
Clonality	Recombinant
NCBI	QHD43423
Dilution Range	Paired with 75G5a, used for specific, high-sensitivity detection of SARS-CoV-2 nucleocapsid protein (NP).
Expiration Date	12 months from date of receipt.