



Product Datasheet

Recombinant SARS-CoV-2 (COVID-19) Nucleocapsid Antibody [6G9] (orb1274615)


Description	Recombinant SARS-CoV-2 (COVID-19) Nucleocapsid Antibody [6G9]
Species/Host	Human
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	full length recombinant Nucleocapsid Protein
Target	N
Preservatives	PBS, pH 7.4
Form/Appearance	Liquid
Concentration	5.0 mg/ml (+/- 10 %)
Storage	24 months from manufacturing at -20 °C; 1 month, 2 to 8 °C under sterile conditions.
Note	For research use only
Application notes	Immunoreactivity: 80-120 % compared to the reference sample in an ELISA test.
Isotype	Human IgG1
Clonality	Recombinant
Uniprot ID	P0DTC9
NCBI	QHD43423
Dilution Range	Immunoreactivity: 80-120 % compared to the reference sample in an ELISA test.
Expiration Date	12 months from date of receipt.



Immobilized
2019-nCoV
Nucleocapsid
Prote...



Loaded 2019-
nCoV
Nucleocapsid
Protein-Hi...



Please note that
pair
recommendations
ar...