

Product Datasheet

SARS-CoV-2 Spike RBD Antibody [2355] (orb1274621)

Description

SARS-CoV-2 Spike RBD Antibody

Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Sequenced from human survivors of COVID-19 (SARS-CoV-2). These are made recombinantly in HEK293 cells. The spike RBD is expressed on the surface of SARS-CoV-2.
Preservatives	This monoclonal antibody is aseptically packaged and formulated in 0.01 M phosphate buffered saline (PBS) pH 7.2 - 7.4, 150 mM NaCl with no carrier protein, potassium, calcium or preservatives added.
Form/Appearance	Liquid
Concentration	> 1.0 mg/ml
Storage	Functional grade preclinical antibodies may be stored sterile as received at 2-8° C for up to one month. For longer term storage, aseptically aliquot in working volumes without diluting and store at -80° C. Avoid Repeated Freeze Thaw Cycles.
Note	For research use only
Application notes	ELISA (Quality Tested by Biorbyt)
Isotype	Human IgG1
Clonality	Recombinant
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558