


Product Datasheet

SARS-CoV-2 Spike NTD Antibody [2146] (HRP) (orb1274624)

Description	SARS-CoV-2 Spike NTD Antibody [2146] (HRP)
Reactivity	Virus
Conjugation	HRP
Tested Applications	ELISA, IHC
Immunogen	Sequenced from human survivors of COVID-19 (SARS-CoV-2). The spike NTD is expressed on the surface of SARS-CoV-2.
Target	S
Preservatives	This HRP-conjugated antibody is formulated in 0.01 M phosphate buffered saline (PBS) pH 7.4, 1% rAlbumin. (Warning: Use of sodium azide as a preservative will inhibit the enzyme activity of horseradish peroxidase)
Form/Appearance	Liquid
Concentration	0.5 mg/ml
Storage	This horseradish peroxidase conjugated monoclonal antibody is stable when stored at 2-8° C. Do not freeze.
Note	For research use only
Application notes	ELISA (Quality Tested by Biobyte); IHC (Other applications as reported in literature)
Isotype	Human IgG1
Clonality	Recombinant
Expiration Date	12 months from date of receipt.



Coating:
Purified
Recombinant
SARS-CoV-
2...