

Product Datasheet

SARS-CoV-2 Spike NTD Antibody [2146] (Biotin) (orb1274625)

Description SARS-CoV-2 Spike NTD Antibody [2146] (Biotin)

Reactivity Virus

Conjugation Biotin

Tested Applications ELISA, IHC

Immunogen Sequenced from human survivors of COVID-19 (SARS-CoV-2). The spike NTD is expressed on the surface of SARS-CoV-2.

Target S

Preservatives This biotinylated monoclonal antibody is formulated in 0.01 M phosphate buffered saline (PBS) pH 7.4, 1% rAlbumin and 0.09% sodium azide as a preservative.

Form/Appearance Liquid

Concentration 0.5 mg/ml

Storage This biotinylated monoclonal antibody is stable when stored at 2-8° C. Do not freeze. Thaw Cycles.


Note For research use only

Application notes ELISA (Quality Tested by Biorbyt); IHC (Other applications as reported in literature)

Isotype Human IgG1

Clonality Recombinant

Expiration Date 12 months from date of receipt.



Coating:
Purified
Recombinant
SARS-CoV-
2...