



## Product Datasheet

### PAK3 antibody (orb127465)

|                            |                                    |
|----------------------------|------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to PAK3 |
| <b>Species/Host</b>        | Rabbit                             |
| <b>Reactivity</b>          | Human, Mouse, Rat                  |
| <b>Conjugation</b>         | Unconjugated                       |
| <b>Tested Applications</b> | WB                                 |
| <b>Immunogen</b>           | Internal                           |
| <b>Concentration</b>       | 1mg/ml                             |
| <b>Note</b>                | For research use only              |
| <b>Clonality</b>           | Polyclonal                         |
| <b>Uniprot ID</b>          | <b>O75914</b>                      |
| <b>Dilution Range</b>      | WB: 1:500-1:3000, ELISA: 1:40000   |
| <b>Expiration Date</b>     | 12 months from date of receipt.    |



Western blot analysis of extracts from r...



Western blot analysis of extracts from H...