



Product Datasheet

KIF5A antibody (orb12749)

Description	Rabbit polyclonal antibody to KIF5A
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human KIF5A(NP_004975.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	129kDa
Uniprot ID	Q12840
Dilution Range	WB: 1:1000-1:3000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse brain...



Western blot analysis of extracts of var...