


Product Datasheet

SLITRK6 antibody (orb12750)

Description	Rabbit polyclonal antibody to SLITRK6
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human SLITRK6(NP_115605.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	120kDa
Uniprot ID	Q9H5Y7
Dilution Range	WB: 1:500-1:2000, IF/ICC: 1:50-1:100, FACS: 1:20-1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...