



## Product Datasheet

### WAVE1 antibody (orb127522)


|                            |                                                                          |
|----------------------------|--------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to WAVE1                                      |
| <b>Species/Host</b>        | Rabbit                                                                   |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                        |
| <b>Conjugation</b>         | Unconjugated                                                             |
| <b>Tested Applications</b> | IF, IHC, WB                                                              |
| <b>Immunogen</b>           | Internal                                                                 |
| <b>Concentration</b>       | 1mg/ml                                                                   |
| <b>Note</b>                | For research use only                                                    |
| <b>Clonality</b>           | Polyclonal                                                               |
| <b>Uniprot ID</b>          | <a href="#">Q92558</a>                                                   |
| <b>Dilution Range</b>      | WB: 1:500-1:3000, IHC-P: 1:50-1:100, IF/ICC: 1:100-1:500, ELISA: 1:40000 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                          |



Western blot analysis of extracts from 293T cells transfected with WAVE1.



Immunohistochemical analysis of formalin-fixed paraffin-embedded (FFPE) sections of human tissue.



Confocal immunofluorescence analysis of cells expressing WAVE1.