

Product Datasheet

PCDHGB5 Antibody (orb1275437)

Description	Rabbit polyclonal antibody to PCDHGB5
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Synthetic peptide of human PCDHGB5
Preservatives	pH 7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Immunogen-specific rabbit IgG
Clonality	Polyclonal
Uniprot ID	Q9Y5G0
NCBI	NP_061748
Dilution Range	IHC: 1/30-1/150; ELISA: 1/5000-1/10000
Expiration Date	12 months from date of receipt.