


Product Datasheet

Histone H2A antibody (orb127582)

Description	Rabbit polyclonal antibody to Histone H2A
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P04908 , Q7L7L0 , P0C0S8 , Q96KK5 , P20671 , Q6FI13
Dilution Range	IHC-P: 1:50-1:100, ELISA: 1:5000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...