

## Product Datasheet

### Integrin beta1 antibody (orb127588)

|                            |                                              |
|----------------------------|----------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to Integrin beta1 |
| <b>Species/Host</b>        | Rabbit                                       |
| <b>Reactivity</b>          | Human, Mouse, Rat                            |
| <b>Conjugation</b>         | Unconjugated                                 |
| <b>Tested Applications</b> | WB                                           |
| <b>Immunogen</b>           | C-terminal                                   |
| <b>Concentration</b>       | 1mg/ml                                       |
| <b>Note</b>                | For research use only                        |
| <b>Clonality</b>           | Polyclonal                                   |
| <b>Uniprot ID</b>          | <b>P05556</b>                                |
| <b>Dilution Range</b>      | WB: 1:500-1:3000, ELISA: 1:1000              |
| <b>Expiration Date</b>     | 12 months from date of receipt.              |

