



Product Datasheet

GOSR1 antibody (orb12760)

Description	Rabbit polyclonal antibody to GOSR1
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human GOSR1(NP_004862.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	29kDa
Uniprot ID	O95249
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U2OS cell...



Western blot analysis of extracts of mou...