

Product Datasheet

SLC4A3 Antibody (orb1276018)

Description	Rabbit polyclonal antibody to SLC4A3
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Synthetic peptide of human SLC4A3
Preservatives	pH 7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Immunogen-specific rabbit IgG
Clonality	Polyclonal
MW	136 kDa
Uniprot ID	P48751
NCBI	NP_005061
Dilution Range	WB: 1/1000-1/5000; IHC: 1/25-1/50; ELISA: 1/5000-1/10000
Expiration Date	12 months from date of receipt.