



Product Datasheet

AARSD1 antibody (orb127616)

Description	Rabbit Polyclonal to AARSD1.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q9BTE6
Dilution Range	WB: 1:500-1:3000, ELISA: 1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from H...



Western blot analysis of extracts from H...