



Product Datasheet

Actinin alpha-2/3 antibody (orb127633)

Description	Rabbit Polyclonal to Actinin alpha-2/3.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	N-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q08043 , P35609
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:100, IF/ICC: 1:100-1:500, ELISA: 1:20000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from H...



Immunohistochemical analysis of formalin...



Confocal immunofluorescence analysis of ...