

Product Datasheet

ABCC8 Antibody (orb1276366)

Description	Rabbit polyclonal antibody to ABCC8
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Synthetic peptide of human ABCC8
Preservatives	pH 7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Immunogen-specific rabbit IgG
Clonality	Polyclonal
Uniprot ID	Q09428
NCBI	NP_000343
Dilution Range	IHC: 1/35-1/150; ELISA: 1/2000-1/5000
Expiration Date	12 months from date of receipt.