


## Product Datasheet

### ALDH1A2 antibody (orb127673)

|                            |                                  |
|----------------------------|----------------------------------|
| <b>Description</b>         | Rabbit Polyclonal to ALDH1A2.    |
| <b>Species/Host</b>        | Rabbit                           |
| <b>Reactivity</b>          | Human, Mouse, Rat                |
| <b>Conjugation</b>         | Unconjugated                     |
| <b>Tested Applications</b> | WB                               |
| <b>Immunogen</b>           | C-terminal                       |
| <b>Concentration</b>       | 1mg/ml                           |
| <b>Note</b>                | For research use only            |
| <b>Clonality</b>           | Polyclonal                       |
| <b>Uniprot ID</b>          | <b>O94788</b>                    |
| <b>Dilution Range</b>      | WB: 1:500-1:3000, ELISA: 1:10000 |
| <b>Expiration Date</b>     | 12 months from date of receipt.  |



Western blot analysis of extracts from K...