

## Product Datasheet

### DUSP2 Antibody (orb1276894)

|                            |                                                                                          |
|----------------------------|------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to DUSP2                                                      |
| <b>Species/Host</b>        | Rabbit                                                                                   |
| <b>Reactivity</b>          | Human, Mouse                                                                             |
| <b>Conjugation</b>         | Unconjugated                                                                             |
| <b>Tested Applications</b> | ELISA, IHC, WB                                                                           |
| <b>Immunogen</b>           | Synthetic peptide of human DUSP2                                                         |
| <b>Preservatives</b>       | pH 7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol                                        |
| <b>Storage</b>             | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |
| <b>Note</b>                | For research use only                                                                    |
| <b>Isotype</b>             | Immunogen-specific rabbit IgG                                                            |
| <b>Clonality</b>           | Polyclonal                                                                               |
| <b>MW</b>                  | 34 kDa                                                                                   |
| <b>Uniprot ID</b>          | <b>Q05923</b>                                                                            |
| <b>NCBI</b>                | <b>NP_004409</b>                                                                         |
| <b>Dilution Range</b>      | WB: 1/500-1/2000; IHC: 1/50-1/200; ELISA: 1/5000-1/10000                                 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                          |