

Product Datasheet

CACNA1G Antibody (orb1277006)

Description	Rabbit polyclonal antibody to CACNA1G
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Synthetic peptide of human CACNA1G
Preservatives	pH 7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Immunogen-specific rabbit IgG
Clonality	Polyclonal
Uniprot ID	O43497
NCBI	NP_061496
Dilution Range	IHC: 1/15-1/50; ELISA: 1/1000-1/2000
Expiration Date	12 months from date of receipt.