

Product Datasheet

ANXA2R Antibody (orb1277032)

Description	Rabbit polyclonal antibody to ANXA2R
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Synthetic peptide of human ANXA2R
Preservatives	pH 7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Immunogen-specific rabbit IgG
Clonality	Polyclonal
Uniprot ID	Q3ZCQ2
NCBI	NP_001014301
Dilution Range	IHC: 1/25-1/100; ELISA: 1/1000-1/2000
Expiration Date	12 months from date of receipt.