

Product Datasheet

AGBL1 Antibody (orb1277450)

Description	Rabbit polyclonal antibody to AGBL1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Synthetic peptide of human AGBL1
Preservatives	pH 7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Immunogen-specific rabbit IgG
Clonality	Polyclonal
Uniprot ID	Q96MI9
NCBI	NP_689549
Dilution Range	IHC: 1/15-1/50; ELISA: 1/1000-1/5000
Expiration Date	12 months from date of receipt.