



Product Datasheet

BHLHB3 antibody (orb127757)

Description	Rabbit Polyclonal to BHLHB3.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q9C0J9
Dilution Range	WB: 1:500-1:3000, ELISA: 1:40000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from R...



Western blot analysis of extracts from 3...