



Product Datasheet

STAT3 antibody (orb12776)


Description	Rabbit polyclonal antibody to STAT3
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Recombinant fusion protein of human STAT3(NP_644805.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	83kDa/87kDa/88kDa
Uniprot ID	P40763
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:100, IF/ICC: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse kidney tissue.



Immunohistochemical staining of mouse spleen tissue.



Immunohistochemical staining of human liver tissue.