



Product Datasheet

Syntaxin 1A antibody (orb12780)


Description	Rabbit polyclonal antibody to Syntaxin 1A
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human STX1A(NP_004594.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	33kDa
Uniprot ID	Q16623
Dilution Range	WB: 1:500-1:2000, IP: 1:50-1:100, FACS: 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of C6 cells ...



Immunohistochemical staining of human gastric tissue ...



Immunohistochemical staining of human liver tissue ...