

Product Datasheet

ARPC5L Antibody (orb1278615)

Description	Rabbit polyclonal antibody to ARPC5L
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Fusion protein of human ARPC5L
Preservatives	pH 7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Immunogen-specific rabbit IgG
Clonality	Polyclonal
Uniprot ID	Q9BPX5
Dilution Range	IHC: 1/30-1/150; ELISA: 1/5000-1/10000
Expiration Date	12 months from date of receipt.