

Product Datasheet

RNF24 Antibody (orb1278618)

Description	Rabbit polyclonal antibody to RNF24
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Fusion protein of human RNF24
Preservatives	pH 7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Immunogen-specific rabbit IgG
Clonality	Polyclonal
MW	17 kDa
Uniprot ID	Q9Y225
Dilution Range	WB: 1/200-1/1000; IHC: 1/20-1/100; ELISA: 1/5000-1/10000
Expiration Date	12 months from date of receipt.