

Product Datasheet

Chk2 (phospho-Thr383) antibody (orb127864)

Description	Rabbit polyclonal antibody to Chk2 (phospho-Thr383).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	O96017
Dilution Range	WB: 1:500-1:3000, ELISA: 1:1000
Expiration Date	12 months from date of receipt.

