


Product Datasheet

APOBEC3C antibody (orb12790)

Description	Rabbit polyclonal antibody to APOBEC3C
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF
Immunogen	Recombinant fusion protein of human APOBEC3C(NP_055323.2).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	22kDa
Uniprot ID	Q9NRW3
Dilution Range	IF/ICC: 1:20 - 50
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U2OS cell...