

Product Datasheet

CIAO2A Antibody (orb1279482)

Description	Rabbit polyclonal antibody to CIAO2A
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Fusion protein of human CIAO2A
Preservatives	pH 7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Immunogen-specific rabbit IgG
Clonality	Polyclonal
MW	18 kDa
Uniprot ID	Q9H5X1
Dilution Range	WB: 1/200-1/1000; ELISA: 1/5000-1/10000
Expiration Date	12 months from date of receipt.